



***ATP Biotech Inc.***

---

**ATP™ Reagent Genomic DNA Kit**  
**Catalog No. AGE100**



# ATP™ Reagent Genomic DNA Kit

Store at room temperature (15-25°C)

## Contents

|  |   |
|--|---|
| Blood (50 ~ 300 $\mu$ l) Protocol .....  | 1 |
| Blood (300 $\mu$ l ~ 3ml) Protocol ..... | 2 |
| Blood (3 ~ 10 ml) Protocol .....         | 3 |
| Culture Cell Protocol .....              | 4 |
| Gram-Positive Bacteria Protocol .....    | 5 |
| Yeast Protocol .....                     | 5 |

1. ATP™ provide RBC Lysis Buffer to remove non-nucleated red blood cells and reduce hemoglobin contamination. But when blood sample is less than 50  $\mu$ l or sample is nucleated red blood cells, we recommend using Cultured Cells protocol to purify genomic DNA.
2. ATP™ provide genomic DNA extraction protocols for smaller sample quantity. If larger sample quantity is required, user can scale up the buffer volume of the protocols proportionally.

## Blood Protocol (50~300 $\mu$ l)

Additionally required : Microcentrifuge tube 、 Isopropanol 、 70 % Ethanol 、 RNase A (10 mg/ml)

### RBC Lysis (Use fresh blood)

1. Collected fresh blood in EDTA-NA2-treated collection tubes (or other anticoagulant).
2. Apply up to 300  $\mu$ l of blood to a 1.5ml microcentrifuge tube. If blood sample is more than 300  $\mu$ l blood (up to 1 ml), apply the blood sample to a sterile 15ml centrifuge tube.
3. Add 3X the sample volume of RBC Lysis Buffer and mix by inversion. Do not vortex.
4. Incubate the tube for 5 min at room temperature.
5. Centrifuge for 2 min at 3000 x g.
6. Remove the supernatant, but retain about 50  $\mu$ l residual Buffer to resuspend the white cell pellet by vortexing.

### Cell Lysis

7. Add 300  $\mu$ l Cell Lysis Buffer to the tube and mix by vortexing.
8. Incubate at 60 °C for 10 minutes until the sample lysate is clear. During incubation, invert the tube every 3 minutes.

Optional Step: RNA degradation (If RNA-free genomic DNA is required, perform this optional step.)

- Add 2  $\mu$ l of RNase A (10 mg/ml, provided by user) to sample lysate and mix by vortexing.
- Incubate at room temperature for 10 minutes.

### Protein Remove

9. Add 100  $\mu$ l of Protein Remove Buffer to the sample lysate and mix immediately by vortexing for 10 seconds.
10. Incubate on ice for 5 minutes.
11. Centrifuge at full speed (13,000 rpm) for 3 minutes.

### DNA Precipitation

12. Transfer the supernatant from Step 11 to a microcentrifuge tube.
13. Add 300  $\mu$ l of Isopropanol and mix well by inverting.
14. Centrifuge at full speed (13,000 rpm) for 3 minutes.
15. Discard the supernatant and add 300  $\mu$ l of 70% ethanol to wash the pellet.
16. Centrifuge at full speed (13,000 rpm) for 1 minutes.
17. Discard the supernatant and air-dry the pellet for 10 minutes.

### DNA Rehydration

18. Add 50-100  $\mu$ l of Water or TE buffer and incubate at 60°C for 30-60 minutes to dissolve DNA pellet. During incubation, tap the bottom of tube to promote DNA rehydration.

### Blood Protocol (300 $\mu$ l ~ 3ml)

Additionally required : 15 ml centrifuge tube 、 Isopropanol 、 70 % Ethanol 、 RNase A (10 mg/ ml)

#### RBC Lysis (Use fresh blood)

1. Collected fresh blood in EDTA-NA2-treated collection tubes (or other anticoagulant).
2. Apply up to blood sample (up to 3ml) to a sterile 15ml centrifuge tube.
3. Add 3X the sample volume of RBC Lysis Buffer and mix by inversion. Do not vortex.
4. Incubate the tube for 5 min at room temperature.
5. Centrifuge at 2000 x g for 5 min.
6. Remove the supernatant, but retain about 300  $\mu$ l residual Buffer to resuspend the white cell pellet by vortexing.

#### Cell Lysis

7. Add 3ml Cell Lysis Buffer to the tube and mix by vortexing.
8. Incubate at 60 °C for 10 minutes until the sample lysate is clear. During incubation, invert the tube every 3 minutes.

Optional Step: RNA degradation (If RNA-free genomic DNA is required, perform this optional step.)  
· Add 10  $\mu$ l of RNase A (10 mg/ml, provided by user) to sample lysate and mix by vortexing.  
· Incubate at room temperature for 10 minutes.

#### Protein Remove

9. Add 1ml of Protein Remove Buffer to the sample lysate and mix immediately by vortexing for 10 seconds.
10. Incubate on ice for 5 minutes.
11. Centrifuge at full speed (13,000 rpm) for 5 minutes.

#### DNA Precipitation

12. Transfer supernatant (about 4 ml) from Step11 to a 15ml centrifuge tube.
13. Add 3ml of Isopropanol and mix well by inverting.
14. Centrifuge at full speed (13,000 rpm) for 5 minutes.
15. Carefully remove the supernatant and add 3ml of 70% ethanol to wash the pellet.
16. Centrifuge at full speed (13,000 rpm) for 1 min.
17. Discard the supernatant and air-dry the pellet for 20 minutes.

#### DNA Rehydration

18. Add 100-300  $\mu$ l of Water or TE buffer and incubate at 60°C for 30-60 minutes to dissolve DNA pellet. During incubation, tap the bottom of tube to promote DNA rehydration.

### Blood Protocol (3ml ~ 10ml)

Additionally required : 50 ml centrifuge tube 、 Isopropanol 、 70 % Ethanol 、 RNase A (10 mg/ ml)

#### RBC Lysis (Use fresh blood)

1. Collected fresh blood in EDTA-NA2-treated collection tubes (or other anticoagulant).
2. Apply up to blood sample (up to 10 ml) to a sterile 50ml centrifuge tube.
3. Add 3X the sample volume of RBC Lysis Buffer and mix by inversion. Do not vortex.
4. Incubate the tube for 5 min at room temperature.
5. Centrifuge at 2000  $\mu$ g for 5 min.
6. Remove the supernatant, but retain about 500  $\mu$ l residual Buffer to resuspend the white cell pellet by vortexing.

#### Cell Lysis

7. Add 10ml Cell Lysis Buffer to the tube and mix by vortexing.
8. Incubate at 60 °C for 10 minutes until the sample lysate is clear. During incubation, invert the tube every 3 minutes.

Optional Step: RNA degradation (If RNA-free genomic DNA is required, perform this optional step.)

- Add 20  $\mu$ l of RNase A (10 mg/ml, provided by user) to sample lysate and mix by vortexing.
- Incubate at room temperature for 10 minutes.

#### Protein Remove

9. Add 3.3 ml of Protein Remove Buffer to the sample lysate and mix immediately by vortexing for 10 seconds.
10. Incubate on ice for 5 minutes.
11. Centrifuge at full speed (13,000 rpm) for 3 minutes.

#### DNA Precipitation

12. Transfer supernatant (about 4 ml) from Step 11 to a 15 ml centrifuge tube.
13. Add 10 ml of Isopropanol and mix well by inverting.
14. Centrifuge at full speed (13,000 rpm) for 5 minutes.
15. Carefully remove the supernatant and add 10ml of 70% ethanol to wash the pellet.
16. Centrifuge at full speed (13,000 rpm) for 1 min.
17. Discard the supernatant and air-dry the pellet for 20 minutes.

#### DNA Rehydration

18. Add 300-600  $\mu$ l of Water or TE buffer and incubate at 60 °C for 30-60 minutes to dissolve DNA pellet. During incubation, tap the bottom of tube to promote DNA rehydration.

## Culture Cells Protocol

Additionally required : Microcentrifuge tube 、 Isopropanol 、 70 % Ethanol 、 RNase A (10 mg/ ml)

### Sample Preparation

#### (Cultured animal cells)

If use adherent cells, trypsinize the cells before harvesting.

- A. Transfer  $10^6$ - $10^7$  of cells to a microcentrifuge tube (not provided) and harvest the cells with centrifugation for 20 seconds at  $6,000 \times g$  (about 8,000 rpm for microcentrifuge).
- B. Remove the supernatant, but retain about 50  $\mu$ l residual buffer to resuspend the white cell pellet by vortexing.

#### (Blood)

For mammalian blood (non-nucleated), the sample volume is up to 50  $\mu$ l.

For nucleated erythrocytes (e.g., bird or fish), the sample volume is up to 10  $\mu$ l.

- A. Collected fresh blood in EDTA-NA2-treated collection tubes (or other anticoagulant).

### Cell Lysis

1. Add 300  $\mu$ l Cell Lysis Buffer to the sample and mix by vortexing.
2. Incubate at 60 °C for 10 minutes until the sample lysate is clear. During incubation, invert the tube every 3 minutes.

Optional Step: RNA degradation (If RNA-free genomic DNA is required, perform this optional step.)

- Add 4  $\mu$ l of RNase A (10 mg/ml, provided by user) to sample lysate and mix by vortexing.
- Incubate at room temperature for 5 minutes.

### Protein Remove

3. Add 100  $\mu$ l of Protein Remove Buffer to the sample lysate and mix immediately by vortexing for 10 seconds.
4. Incubate on ice for 5 minutes.
5. Centrifuge at full speed (13,000 rpm) for 3 minutes.

### DNA Precipitation

6. Transfer the supernatant from Step5 to a microcentrifuge tube.
7. Add 300  $\mu$ l of Isopropanol and mix well by inverting.
8. Centrifuge at full speed (13,000 rpm) for 3 minutes.
9. Discard the supernatant and add 300  $\mu$ l of 70% ethanol to wash the pellet.
10. Centrifuge at full speed (13,000 rpm) for 3 minutes.
11. Discard the supernatant and air-dry the pellet for 10 minutes.

### DNA Rehydration

12. Add 50-100  $\mu$ l of Water or TE buffer and incubate at 60 °C for 30-60 minutes to dissolve DNA pellet. During incubation, tap the bottom of tube to promote DNA rehydration.

## Gram-positive Bacterial Protocol

If sample is gram-negative bacteria, use Cultured Cells Protocol.

Additionally required : Lysozyme Buffer (20 mg/ml lysozyme ; 20 mM Tris-HCl ; 2 mM EDTA ; 1% Triton X-100 , pH 8.0) , Prepare the lysozyme buffer just before use

### Cell Harvest / Prelysis (For Gram-positive bacteria)

1. Transfer bacterial culture ( $< 10^9$ ) to a microcentrifuge tube (not provided).
2. Centrifuge for 1 min at full speed (13,000 rpm) in a microcentrifuge and discard the supernatant.
3. Add 100  $\mu$ l of Lysozyme Buffer to the tube and resuspend the cell pellet by vortexing or pipetting.
4. Incubate at room temperature for 20 minutes. During incubation, invert the tube every 2-3min.
5. Proceed with Step1 Cell Lysis of Cultured Cells protocol.

## Yeast Protocol

Additionally required : Sorbitol buffer (1.2M sorbitol ; 10mM  $\text{CaCl}_2$  ; 0.1M Tris/Cl , pH 7.5 ; 35mM mercaptoethanol)

Lyticase or zymolase

### Cell Harvest / Prelysis (For Gram-positive bacteria)

1. Harvest yeast cells (up to  $5 \times 10^7$ ) by centrifugation for 10 min at 5,000g.
2. Resuspend the pellet in 600  $\mu$ l sorbitol buffer.
3. Add 200U of lyticase or zymolase. Incubate at 30 °C for 30 min.
4. Centrifuge the mixture for 10 min at 2,000g to harvest Spheroplast.
5. Proceed with Step1 Cell Lysis of Cultured Cells protocol.

