



ATP Biotech Inc.

ATP™ Viral Nucleic Acid Extraction Kit
Catalog No. AVR050/AVR100



ATP™ Viral Nucleic Acid Extraction Kit

Store at room temperature (15-25°C)

Contents

ATP™ Viral Nucleic Acid Extraction Kit

Introduction	1
Kit Contents	1
Cautions	2
Use Limitation	2
Recommend	2
Equipments and Reagents are provided by User	2
Viral Nucleic Acid Extraction Protocol	3

Introduction

Format : Spin column
 Sample : 200 µl of serum, plasma, body fluid and cell culture supernatant
 Operation : Centrifuge / vacuum manifold
 Operation time : 40 minutes
 Elution volume : 50 µl
 Application : RT-PCR, PCR, Real-Time PCR, Real-Time RT-PCR, Automated Fluorescent DNA Sequencing

ATP™ Viral Nucleic Acid Extraction Kit is specially designed for high-throughput purification of viral RNA or DNA from cell-free samples as serum, plasma, body fluids, and the supernatant of viral infected cell culture. With the extraction method, DNA/RNA viruses are lysed quickly and efficiently by lysis buffer which is a highly concentrated solution of chaotropic salt. Nucleic acid in chaotropic salt and methanol are bond to the glass-fiber matrix of viral DNA/RNA Binding Column. Contaminations like salts, metabolites and soluble macromolecular components are removed in Washing Steps. The nucleic acids can be eluted by low-salt buffer or water and are ready-to-use in subsequent reactions. The detection limit for certain viruses depends on the sensitivity of individual PCR or RT-PCR assay. This protocol is recommended for parallel purification of viral RNA including HCV, HIV, HTLV and viral DNA including HBV and CMV.

Kit Contents : Cat.No. / Kit Contents

AVR050 (50 preps/kit)	AVR100 (100 preps/kit)
VB Lysis Buffer : 30 ml	VB Lysis Buffer : 60 ml
AD Buffer (concentrated)* : 4 ml	AD Buffer (concentrated)*** : 8 ml
W1 Buffer : 30 ml	W1 Buffer : 50 ml
Wash Buffer (concentrated)** : 12.5 ml	Wash Buffer (concentrated)**** : 25 ml
RNase-free water : 6 ml	RNase-free water : 6 ml
VB Columns : 50 pcs (yellow/white filter)	VB Columns : 100 pcs (yellow/white filter)
2ml Collection Tubes : 100 pcs	2ml Collection Tubes : 200 pcs

* Add 30 ml ethanol (96-100%) to AD Buffer prior to initial use.(see bottle label for volume)
 ** Add 50 ml ethanol (96-100%) to Wash Buffer prior to initial use.
 *** Add 60 ml ethanol (96-100%) to AD Buffer prior to initial use. (see bottle label for volume)
 **** Add 100 ml ethanol (96-100%) to Wash Buffer prior to initial use.

Caution : VB Lysis Buffer contain guanidine hydrochloride which is a harmful and irritant agent. During operation, always wear a lab coat, disposable gloves, and protective goggles.

Product Intended Use : Research / Clinical Application

ATP™ Viral Nucleic Acid Extraction Kit is a research purpose device. ATP™ Biotech Inc. has not validated in clinical application for any organism or association, and therefore offer no specific claims for uses in diagnostics or prognostics. This device can serve as a means for molecular assays in clinical diagnostics laboratory systems after the laboratory has certified their systems according to the CLIA¹88 regulation in the USA or local equivalents in other countries. Exercise all necessary care and attention when handling this product.

Notice :

ATP™ Viral Nucleic Acid Extraction Kit could extract both DNA and RNA from samples containing DNA and RNA based viruses. When user do RNA virus extraction, it is essential to use **cell free** samples to reduce DNA contamination.

Equipments and Reagents are provided by User

- 1.5 ml (RNase-free) microcentrifuge tubes
- Microcentrifuge with rotor for 2 ml tubes
- Ethanol (96-100%)
- PBS (phosphate-buffered saline)

Viral Nucleic Acid Extraction Kit Protocol

- AVR050 : Add 30 ml ethanol (96-100%) to AD Buffer prior to the initial use. Add 50 ml ethanol (96-100%) to Wash Buffer prior to the initial use.
- AVR100 : Add 60 ml ethanol (96-100%) to AD Buffer prior to the initial use. Add 100 ml ethanol (96-100%) to Wash Buffer prior to the initial use.
- Additionally required : PBS \ Ethanol(96-100%) \ 1.5 ml microcentrifuge tube (RNase-free)

Lysis

1. Transfer 200 μ l of sample (serum, plasma, body fluids and the supernatant of viral infected cell culture) into a microcentrifuge tube (provided by user). If prepared sample is less than 200 μ l, adjust sample volume to 200 μ l with PBS (provided by user).
2. Add 400 μ l of VB Lysis Buffer to the sample, mix by vortexing.
3. Incubate at room temperature for 10 minutes.

Nucleic Acid Binding

4. Place a VB Column in a 2 ml Collection Tube.
5. Add 450 μ l of AD Buffer (ethanol added) to the sample lysate and mix immediately by vortexing.
6. Apply 600 μ l of lysate mixture from previous step into the VB column.
7. Centrifuge at 13,000 rpm for 1 minute.
8. Discard the flow-through waste and apply the rest of lysate mixture into the same Column.
9. Centrifuge at 13,000 rpm for 1 minute.
10. Discard the Collection tube containing the flow-through waste and transfer the VB column in a new 2 ml Collection tube.


Washing

11. Add 400 μ l of W1 Buffer into the VB Column.
12. Centrifuge at 10,000 xg (13,000 rpm) for 30 seconds.
13. Discard the flow-through and place the VB Column back in Collection Tube.
14. Add 600 μ l of Wash Buffer (ethanol added) into the VB Column
15. Centrifuge at 10,000 xg (13,000 rpm) for 30 seconds.
16. Discard the flow-through and place the VB Column back in the Collection Tube.
17. Centrifuge again for 3 minutes at full speed (13,000 rpm) to dry the column matrix.

Nucleic Acid Elution

18. Place dried VB column in a clean microcentrifuge tube (RNase-free, provided by user).
19. Apply 50 μ l of RNase-free water onto the center of the VB column matrix, and close the cap of VB Column.
20. Stand for 2 minutes until water is absorbed by the matrix.
21. Centrifuge at full speed for 1 minute to elute purified nucleic acid.

ATP™ Viral Nucleic Acid Extraction Kit



Note

Catalog

Product Name	Package	Cat. No.
ATP™ Plasmid Mini Kit	100/300 prep	APD100/APD300
ATP™ Plasmid Midi Kit (Ultra Pure)	25 prep	API25
ATP™ Plasmid Maxi Kit (Ultra Pure)	10/25 prep	APM10/APM25
ATP™ 96-Well Plasmid Mini Kit	4/10 plates	APD9604/APD9610
ATP™ Plasmid Mini Binding Column	50 pcs	PDC50
ATP™ Plasmid Midi Resin Column	10 pcs	PIC10
ATP™ Plasmid Maxi Resin Column	10 pcs	PMC10
ATP™ Gel/PCR DNA Fragments Extraction Kit	100/300 prep	ADF100/ADF300
ATP™ 96-Well Gel/PCR DAN Extraction Kit	4/10 plates	ADF9604/ADF9610
ATP™ 96-Well SEQ Dye Clean Up Kit	4/10 plates	ADC9604/ADC9610
ATP™ Fragment DNA Binding Column	50 pcs	DFC50
ATP™ Genomic DNA Mini Kit (Blood/Culture Cell/Bacteria)	100/300 preps	AGB100/AGB300
ATP™ Genomic DNA Mini Kit (Tissue)	50/300 preps	AGT050/AGT300
ATP™ Genomic DNA Mini Kit (Plant)	100 preps	AGP100
ATP™ Genomic DNA Maxi Kit (Fresh Blood)	25 prep	AGBM25
ATP™ Genomic DNA Maxi Kit (Frozen Blood)	25 prep	AGDM25
ATP™ Plant Genomic DNA Maxi Kit	25 prep	AGPM25
ATP™ 96-Well Genomic DNA Kit	4/10 plates	AGB9604/AGB9610
ATP™ Reagent Genomic DNA Kit	For 100ml blood	AGE100
ATP™ Genomic DNA Binding Column	50 pcs	GDC50
ATP™ RNA Mini Kit (Blood/Culture Cell/Bacteria)	50 prep	ARB050
ATP™ RNA Mini Kit (Tissue)	50 prep	ART050
ATP™ RNA Mini Kit (Plant)	50 prep	ARP050
ATP™ Viral Nucleic Acid Mini Kit	50 prep	AVR050
ATP™ 96-Well Viral Nucleic Acid Kit	4/10 plates	AVR9604/AVR9610
ATP™ RNA Maxi Kit	10 prep	ARTM10
ATP™ Plant RNA Maxi Kit	10 prep	ARPM10
ATP™ RNA Binding Column	50 pcs	RBC50
Proteinase K	11 mg/kit	APK000011
RNase A	0.2 ml (50 mg/ml)	ARA500200
RNase A	1.5 ml (50 mg/ml)	ARA501500

